MACHEREY-NAGEL

Urine analysis



URYXXON® 500

- Quick Reads 400 strips per hour
- Easy Touch screen operation
- Hygienic Automatic strip detection



High performance urine analysis

The URYXXON® 500 is an automatic reader for URYXXON® Stick 10 urine test strips. With a capacity of 400 strips per hour, it is ideal for use in hospitals and practices. The "easy-to-use-features" allow hygienic operation with hardly any training.

The URYXXON® 500 provides dependable urine status results to detect early stages of many diseases such as urinary tract infections, kidney diseases or diabetes. The device eliminates the subjectivity of visual strip evaluation and minimizes risks associated with manual result transcription.

The URYXXON® 500 makes urine analysis quicker, easier and more reliable.



Technical specifications

Instrument memory	500 patient test results including name or patient ID						
Sediment	30 user defined sediment parameters						
Interface	User:	Touchscreen display, alphanumeric input, password protection					
	Computer:	1 x USB (type B) and 2 x RS232 interfaces for connection to PC					
		PS/2 interface for keyboard / barcode reader					
Power requirements	110–240 V AC, automatic						
Dimensions/Weight	240 x 280 mm x 150 mm (D x W x H)						
Operation	e: 5–40 °C (41–104 °F)						
	Humidity range:	20-80 % relative humidity, non condensing					
	Calibration:	automatic, self calibrating					
Capacity	400 strips per hour						
CE	CE labeling according to IVDD 98/79/EC						

URYXXON® 500

Field of application

The URYXXON® 500 is intended for the automatic evaluation of URYXXON® Stick 10 urine test strips. It provides reliable results for the semi-quantitative detection of the following parameters in urine: blood, urobilinogen, bilirubin, protein, nitrite, ketones, glucose, pH, specific gravity and leukocytes. URYXXON® 500 and URYXXON® Stick 10 are CE marked according to the IVD-directive 98/79/EC.

Improvements for your lab

Easy and intuitive

Operating the URYXXON® 500 is easy. Access all functions using the easy-to-follow menu on the touch screen. It nearly eliminates the need for training and makes urine analysis more reliable.



Informative and unique

Add microscopic evidence to the urine chemistry results using the sediment options. Powerful user defined filter functions only admit those samples to sediment analysis that need further diagnostic.



Clean and hygienic

Start the measurement without touching the instrument. Simply place the strip on the transport tray, the detector will recognize it and automatically start the measurement.



Fast and efficient

Once turned on, the URYXXON® 500 is immediately ready for use. It reads up to 400 strips per hour ensuring fast urine chemistry results.



Accurate and convenient

Use a barcode reader to enter patient IDs. Set up lists of samples with ease and convenience and quickly add urine chemistry results to the patient file in your laboratory information system.



Reliable and safe

Use the QC menu to maintain clinical reporting standards with ease and convenience. Medi-Test Control solutions are ideal to ensure correct results at any time.



Small and practical

Free up space with the smallest instrument in its class. With a size of only 28 x 24 cm (11 x 9,4 in) the URYXXON® 500 fits onto every bench. The clear design allows access to all parts and makes cleaning easy.



Medi-Test ordering information

strips s per p pack p	930 20	Glucose 1)	Glucose	Ketones	Ascorbic acid	> .⊑		tes					gen	
	930 20			_	Aso	Protein/ Albumin	Blood	Leukozytes	Nitrite	Density	pH value	Bilirubin	Urobilinogen	Oreatinine
03025 0		4)												
30023	930 27	Glucose / Ketone 1)												
93004 9		Protein 2 1)												
93005 9	3028	Ketone 1)												
93006 9	3029	Nitrite 1)												
93012 -	-	Urbi												
93015 9	93037	Combi 2 1)												
93007 9	93030	Combi 3 A®												
93009 9	93032	Combi 5												
93035 9	93036	Combi 5 N [®]												
93055 -	-	Combi 5 S												
93018 9	93078	Combi 6												
93013 9	93034	Combi 6 A												
93010 9	93022	Combi 7												
- 9	3021	Combi 8 L												
93011 9	93023	Combi 9®												
- 9	93056	Combi 10 [®]												
93079 9	93058	Combi 10 [®] L												
- 9	93060	Combi 11												
93077 9	93067	Combi 10 [®] SGL												
- 9	93068	URYXXON® Stick 10 ^{2) 3)}												
930874 -	-	Microalbumin 3) 5)												

Further Medi-Test urine test strips

Test strips	for veterinary	applications						
_	930870	Combi 10 [®] VET 4)						

Test strips for detection of urine adulteration

93019 Medi-Test Adulteration Stick⁴⁾ for the detection of urine adulteration prior to drugs-of-abuse tests with test fields for Creatinine, Glutaraldehyde, Nitrite, pH, Oxidants and Density

All products (except 4) are CE-marked according to the directive 98/79/EC

¹⁾ suitable for patient self-testing 2 for evaluation with reflectometer URYXXON® 500 3 for evaluation with reflectometer URYXXON® Relax 4 not an IVD product (no CE-mark) 3 pack of 24 test strips

Instrument	Instruments for evaluation of urine test strips URYXXON® Stick 10					
93088	URYXXON® Relax, suitable for the evaluation of MT URYXXON® Stick 10					
930080	URYXXON® 500, suitable for the evaluation of MT URYXXON® Stick 10					
Accessorie	s and spare parts					
93038	Medi-Test Control, solution to check Medi-Test urine test strips and URYXXON® instruments					
93071	Printer paper for URYXXON® 300 and URYXXON® 500, pack of 5 rolls					
93074	Barcode scanner for URYXXON® instruments					

www.mn-net.com

MACHEREY-NAGEL



 $\textbf{MACHEREY-NAGEL GmbH \& Co. KG} \cdot \textbf{Neumann-Neander-Str. 6-8} \cdot 52355 \ \textbf{D\"{u}ren} \cdot \textbf{Deutschland}$ CH:

Tel.: +33 388 68 22 68 +49 24 21 969-0 Tel.: +41 62 388 55 00 Tel.: +33 388 51 76 88 +49 24 21 969-199 +41 62 388 55 05 Fax: Fax: Fax: E-Mail: sales-ch@mn-net.com E-Mail: sales-fr@mn-net.com E-Mail: info@mn-net.com



© Sven Hoppe (S. 3) – stock.adobe.com KATEN100015 Medi-Test URYXXON® 500 en
2 $/\,0/0\,/\,11.2019$ DD \cdot Printed in Germany