

TECHNICAL DATA SHEET

Cell and Tissue Culture Dishes slides

codes D200035, D200060, D200100, D200150

code vented	diameter	surface area cm ²	sterile	case quantity	case weight (kg)	case volume (m ³)
D200035	32.8	8.5	STERILE R	96 x 10	5.00	0.028
D200060	52	21,2	STERILE R	60 x 10	10.50	0.050
D200100	87.8	60.8	STERILE R	30 x 10	10.08	0.076
D200150	135.5	143	STERILE R	120 x 1	8.20	0.082



USE

Cell and tissue culture dishes are ideal containers for cell growth and yields point to a small and medium volume of input, and are also useful in sample separation, pretreatment, storage, etc.

Treated surface.

The uniform thickness of the flat bottom wall guarantees a bottom without distortion.

Sterilized by radiation.

Not pyrogenic.

WRITTEN	VERIFIED	AUTHORISED
Cristina Rojas Marketing Department	Anna Mir Technical, Development, Quality and Environment Director	Mónica Torras Sales & Marketing Director